Troubleshooting

If a malfunction occurs, first check the points below before you send the blood pressure monitor in to be repaired. If error codes are shown on the display, proceed as follows:

Switch the measurement computer to sleep mode by pushing the SET button, the user 1 button or the user 2 button (see fig. A). Observe the solutions provided in the following table.

Problem

The error code "E1" and the pulse indicator appear on the display.

Causes and solutions

The pulse is not being measured correctly.

-- Fit the cuff correctly on your arm.

-- Repeat the measurement.

Problem

The error code "E2" appears on the display.

Causes and solutions

You spoke or moved during the measurement.

-- Repeat the measurement without speaking or moving.

Problem

The error code "E3" appears on the display.

Causes and solutions

The cuff is not on your upper arm.

-- Fit the cuff correctly on your arm.

-- Repeat the measurement.

Problem

The error code "E4" appears on the display.

Causes and solutions

A measurement error has occurred.

-- Repeat the measurement.

Problem

The error code "E5" appears on the display.

Causes and solutions

A problem has occurred when inflating the cuff (inflation pressure > 300 mmHg).

-- Make sure that the cuff hose is not kinked or clamped.

-- Fit the cuff correctly on your arm.

-- Repeat the measurement.

Problem

The display is blank, or the battery indicator and the error code "E6" appear on the display.

Causes and solutions

The batteries were not inserted correctly.

-- Remove the batteries and put them back in with the correct polarity.

Batteries are empty.

-- Replace the empty batteries with new ones.

Problem

The error code "E7" and the Bluetooth® indicator appear on the display.

Causes and solutions

Bluetooth® transfer is not possible.

Once you have prepared everything for automatic transmission, repeat the measurement or transfer the measurements manually.

Bluetooth® transfer is not possible.

The Bluetooth® function of your smartphone or tablet is switched off.

-- Switch the Bluetooth® function of your smartphone or tablet on.

The Bluetooth® transmitter of the blood pressure monitor is not working properly.

-- Take the batteries out of the battery compartment.

-- Put the batteries back in.

-- Switch the Bluetooth® function of your smartphone or tablet off and back on.

Your smartphone or tablet is not working properly.

-- Switch off your smartphone or tablet off completely and then back on.

Problem

The error code "E7" and the Bluetooth® indicator appear on the display.

Causes and solutions

The app was not successfully installed on your smartphone or tablet.

-- Reinstall the app and log in with your user data.

Your smartphone or tablet does not support Bluetooth® Smart.

-- Exchange your smartphone or tablet with a model that supports Bluetooth® Smart.

Problem

Not all measurements saved are transmitted.

Causes and solutions

You performed more than 60 measurements without connecting to the app.

The blood pressure monitor can save a total of 60 measurements per user. Afterwards, the first measurement will be overwritten by the last one.

-- Connect the blood pressure monitor with the app on a regular basis to prevent losing measurements.

Use the FAQs (frequently asked questions) on our website as an additional help for your blood pressure monitor: <http://www.krippel-watches.com>

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